


Collect stool (poop) for ova and parasite (O&P) kit

1 Label containers



(your name, date of birth, collection date & time)


2 Wash hands



3 Follow steps 3-5 for bulk stool




4 Unscrew cap. Use spoon on lid to scoop stool from bulk container into O&P container.




(Fill to red line on tube label)


5 Close lid tightly




6 Repeat with second tube




7 Put container in biohazard bag



Keep cool, not frozen



8 Throw out extra stool and supplies



9 Wash hands

10 Call _____ at (____)____-____ to return stool for testing

Questions/Comments?
www.COFoodSafety.org