



# Collect stool (poop) using a swab

**1** Label container



(your name, date & time of sample)

**2** Wash hands



**6** Insert swab into tube



(make sure it touches gel in bottom)

Close lid tightly

**7** Put container in biohazard bag



Keep cool, not frozen


**8** Throw out extra stool and supplies



**9** Wash hands

**10** Call \_\_\_\_\_ at  
(\_\_\_\_)\_\_\_\_-\_\_\_\_ to return  
stool for testing

**4** Take inside cap off and throw out



**5** Cover this part (end of q-tip) with stool

